

# MiniPro

VERDOT



**Configurable for pilot or small scale production**

**Designed for easy cleaning and maintenance**

**Configured for simultaneous in-line dilution and gradient operations for high performance and reduced cost**

**User-friendly software meets U.S. FDA 21 CFR Part 11 requirements**

**Includes complete documentation for regulatory submission**

## Pilot scale or Small scale GMP manufacturing from a benchtop

The MiniPro system by VERDOT combines versatility, precision and rock-solid performance in a benchtop chromatography system.

The system is configurable as:

- Pilot scale with flow range 1-60 LPH
- Small scale GMP manufacturing processes with flow range 5-120 LPH



**MiniPro chromatography system in Process configuration**

### Complete System

The MiniPro system by VERDOT Ips<sup>2</sup> includes all hardware and software necessary in a stand-alone system. A standard laptop PC presents the complete user interface in a graphical, easy to learn system. SCADA software is used to run the chromatography process, to collect batch data (audit trail) and to communicate with the PLC.

The stainless steel cabinet with closed design allows easy cleaning in the process environment.

### Modular Concept

The flexibility of VERDOT Ips<sup>2</sup> design offers the possibility to adapt the MiniPro chromatography system to multiple processes. In addition to the standard configuration, the following options are available (see Fig.1):

- 4 additional outlet valves
- Pressure sensor (after column)
- Air sensor (after bubble trap)
- Additional pH conductivity and temperature sensors (before column)
- Additional refractometer can be installed as an option
- Integration possible with fraction collector

Each option includes the required combination of hardware and software modules.

All standard and optional components have been validated on the VERDOT Ips<sup>2</sup> test platform.

### Wide flow range

The flexible design of the equipment allows to work with multiple flow rates and column diameters, with pilot or process configuration, as shown in Table 1.


With pilot configuration, the pumps are volumetric and can achieve high precision flowrate with no need of flowmeter control.

With process configuration, the membrane technology pumps can achieve high precision flowrate with flowmeters control.

**Table 1. MiniPro System vs Standard PilMod columns by VERDOT**


Column diameter, mm	Linear speed, cm/hr	
	With pilot configuration	With process configuration
100	13-764	64-1528
140	6-390	32-780
200	3-191	16-382

**Fig. 1. PID of the MiniPro system**  
Available options are marked in green



\* Flowmeters are only implemented with process configuration. The volumetric pumps of the pilot configuration don't require flowmeters

## System controller screen



## User-friendly software

The system controller employs a simple, user-friendly interface for data input and programming commands.

The process system is password protected (four customizable access levels), and all events and actions are recorded in accordance with cGMP guidelines.

The software allows the system to be operated in manual or automatic mode. The automatic mode includes:

- Multiple steps
- Configurable fluid paths
- Control of valves
- Multiple end-of-step conditions
- Pausing and halting alarms
- Interactive pausing steps
- Compliant with 21 CFR Part 11

The standard equipment includes full trend review, printing and data export from the system.

Pre-programmed routine for testing and calculating Asymmetry and HETP values is included.

Key data can be monitored remotely through OPC UA.

## Specifications

Nominal Flow rate

Pilot: 1.5-60 L/h  
25-1000 mL/min  
Process: 5-120 L/h  
83.3-2000 mL/min

Operating pressure

Up to 6 bar at column inlet

Flowmeter accuracy

± 2% of the measured value

Dilution accuracy

± 3% of the measured value

Tubing after pumps

ID 4.8 mm

Temperature range

2°C-30°C (60°C for CIP)

Pumps

Pilot: Piston pumps

Process: Membrane pumps

UV wavelength range

190-500 nm

Conductivity range

0-500 mS/cm 0-14

pH range

0-14

Pressure sensor

-1 to 9 bar

Bubble trap

50mL (pilot) / 200mL (process)  
with:

- 2 level capacitive sensors
- drain & air injection valves

Valves

Pneumatic membrane valves

Static Mixer

Stainless steel

Air sensor

Optical fiber technology

Fittings

Sanitary fittings

Material and certificates

Plastic parts: PEEK, PFA, FEP

Elastomers: EPDM, PTFE

Stainless steel: 316 L

Bubble trap: Glass

Gasket: FDA compliant, USP Class VI

Surface finish

Ra < 0.4 µm

Degree protection

IP54

Weight

135 kg

Dimensions (H x W x D)

808 x 646 x 1060 mm

Power supply

220V-1p-50/60Hz  
(or 115V-1p-50/60Hz)

Air supply

6 bar

Control

PLC: OMRON

SCADA: iFIX 5.8 on Laptop

Support

4 feet, installation on a lab table